



Interlagos Retractor System

The Altus Spine **Interlagos** Retractor System is a pedicle based, patient mounted, posterior lumbar access system designed to support MIS TLIF procedures. Patient-mounted blades with pedicle anchors that allow for intraoperative distraction and compression. With the unique toe and telescoping capability of the medial blade, **Interlagos** provides a highly versatile access solution for 1- and 2-level procedures.



Refer to Product Insert (PI-003) for complete system description, warnings, and precautions.



©2020 Altus Spine
U.S. and International Patents Pending
BR-1801 Rev B

1340 Enterprise Drive
West Chester, PA, 19380
Office: (610) 355-4156
Fax: (610) 300-3049
info@altus-spine.com

www.altus-spine.com