




Monaco Pedicle Screw System

The Altus Spine **Monaco** Pedicle Screw System has a low profile which accepts 5.5mm and 6.0mm CoCr and titanium rods. The bone screw thread design optimizes thread purchase and pull out strength, and has a friction fit with the screw head so it remains in the appropriate position for rod implantation. The **Monaco** Pedicle Screw System is designed to maximize locking power of the screw and rod construct.



Refer to Product Insert (PI-010) for complete system description, warnings, and precautions.



©2020 Altus Spine

BR-1401 Rev B

1340 Enterprise Drive
West Chester, PA, 19380
Office: (610) 355-4156
Fax: (610) 300-3049
info@altus-spine.com

www.altus-spine.com