



# Monaco HA Pedicle Screw System

The Altus Spine **Monaco HA** Pedicle Screw System has a low profile which accepts 5.5mm and 6.0mm CoCr and titanium rods. The **Monaco HA** screws have hydroxyapatite (HA) coating over the entire thread length. The HA coating largely consists of calcium and phosphorus, and is intended to enhance the fixation between the pedicle screw and the surrounding bone. In addition, the bone screw thread design optimizes thread purchase and pull out strength, and has a friction fit with the screw head so it remains in the appropriate position for rod implantation.



Refer to Product Insert (PI-010) for complete system description, warnings, and precautions.



©2021 Altus Spine

BR-1402 Rev A

1340 Enterprise Drive  
West Chester, PA, 19380  
Office: (610) 355-4156  
Fax: (610) 300-3049  
info@altus-spine.com

[www.altus-spine.com](http://www.altus-spine.com)