




Sochi OCT System

The Altus Spine **Sochi** OCT System provides a seamless transition from occiput to sacrum. The **Sochi** OCT System is comprised of polyaxial screws, hooks, rods, locking screw assemblies and connectors. The system components may be rigidly locked together in a variety of configurations to promote fusion for a wide variety of patient anatomies.



Refer to Product Insert (PI-013) for complete system description, warnings, and precautions.



©2020 Altus Spine

BR-1000 Rev A

1340 Enterprise Drive
West Chester, PA, 19380
Office: (610) 355-4156
Fax: (610) 300-3049
info@altus-spine.com

www.altus-spine.com