



Valencia MIS Pedicle Screw System

The Altus Spine **Valencia MIS** Pedicle Screw System is a simple, effective pedicle fixation system that incorporates reliable rod delivery with dependable rod reduction. In less invasive procedures, the integrated extension tabs provide additional control and guidance during construct placement and locking.

Removable constraint collar provides additional control

Visual guidance for rod delivery and construct locking

Reliable, over-the-skin rod measurement

30mm integrated rod reduction capability

Robust hybrid square thread minimizes cross threading during rod reduction

Cannulated progressive multi-diameter thread form screws

Refer to Product Insert (PI-001) for complete system description, warnings, and precautions.

