



AltusSpine | Silverstone-P IBF System


The Altus Spine **Silverstone-P** IBF System features a state-of-the-art, engineered, self-distracting bulleted nose for optimal insertion. The bi-convex, anatomically-shaped endplate interface maximizes stability. Implants are available in precision machined PEEK and titanium that incorporates proprietary surface texturing intended to promote bone growth.



Titanium Implants
22mm, 26mm and 30mm lengths



PEEK Implants
22mm, 26mm and 30mm lengths



Implants available
in 6mm – 16mm
heights in 1mm
increments

Refer to Product Insert (PI-009) for complete system description, warnings, and precautions.

