





Imola

Lateral IBF System



The Altus Spine **Imola** Lateral IBF System includes a unique, proprietary exposure method that allows for direct access and visualization of the surgical site. The innovative retractor blades convert directly from hand-held to table-mount stabilized via in situ retractor frame attachment without compromising the established access exposure.



Refer to Product Insert (PI-009) for complete system description, warnings, and precautions.

*Covered by U.S. Patent No.: 9,408,598. Other U.S. and foreign patents additionally pending

© 2017 Altus Spine
U.S. and International Patents Pending

BR-1701 Rev B

**AltusSpine**
www.altus-spine.com

1340 Enterprise Drive
West Chester, PA, 19380
Office: (610) 355-4156
Fax: (610) 300-3049
info@altus-spine.com