


AltusSpine | Brickyard ALIF System

The Altus Spine **Brickyard** ALIF System is designed as a simple, yet complete, anterior lumbar IBF system. Its aggressive teeth prevent implant migration and maintains optimal placement. Implants are available in precision machined PEEK and titanium that incorporates proprietary surface texturing intended to promote bone growth.





Insertor

Large axial graft window

Available footprints:
32mm x 24mm
38mm x 27mm
42mm x 30mm

Available heights:
9mm–19mm
in 2mm increments

Lordotic angles:
8°, 12°, 16°, 20°
and 25°



Insertor

Surface treatment technology

Refer to Product Insert (PI-009) for complete system description, warnings, and precautions.



©2025 Altus Spine

BR-1703 Rev C

1340 Enterprise Drive
West Chester, PA, 19380
Office: (610) 355-4156
Fax: (610) 300-3049
info@altus-spine.com

www.altus-spine.com